

# RAB2B antibody



Catalog Number: 114428

## Product name

RAB2B antibody

## Immunogen

[Human RAB2B Recombinant protein \(GST tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

RAB2B Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HeLa cells. Observed molecular weight by Western-blot: 24kd

## Preparation

This antibody was obtained by immunization of RAB2B recombinant protein (Accession Number: NM\_032846). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114428(RAB2B antibody) at dilution of 1:800