

Anti-ANGPTL1 antibody



Catalog Number: 176428

Product name

WB, IHC-P

Anti-ANGPTL1 antibody

Dilutions

Specificity

WB: 1:500-1:2000

Human, Mouse, Rat

IHC-P: 1:100-1:300

Antibody description

Rabbit polyclonal antibody to ANGPTL1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human angptl1. aa range:11-60

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

Validations




Fig1:: Western Blot analysis of K562 cells using Angptl1 Polyclonal Antibody.. Secondary antibody [catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded mouse-heart, antibody was diluted at 1:100