

Mouse IFNA2 recombinant protein



Catalog Number: 500828

General Information

Gene Name Synonym

MCG123693

Protein Construction

A DNA sequence encoding the mouse Ifna2 (NP_034633.2) (Cys24-Glu190) was expressed.

Organism

Mouse

Expression Host

Yeast

QC Testing

Activity

Measured in antiviral assays using L929 cells infected with vesicular stomatitis virus. The ED₅₀ for this effect is typically 9-36 pg/mL.

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Cys 24

Molecular Mass

The recombinant mouse Ifna2 consists 167 amino acids and predicts a molecular mass of 19.4 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse IFNA2 / Interferon alpha 2 Protein