

# GAL3ST1 antibody



Catalog Number: 110828

## Product name

GAL3ST1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

GAL3ST1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HEK-293 cells, SGC-7901 cells. Observed molecular weight by Western-blot: 48kd

## Preparation

This antibody was obtained by immunization of GAL3ST1 recombinant protein (Accession Number: NM\_001318103). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110828(GAL3ST1 antibody) at dilution of 1:500