

Canine CD137 / 4-1BB (His Tag) recombinant protein



Catalog Number: 504828

General Information

Protein Construction

A DNA sequence encoding the canine TNFRSF9 (XP_850336.1) (Met1-Ser185) was expressed with a C-terminal polyhistidine tag.

Organism

Canine

Expression Host

Human Cells

QC Testing

Activity

Measured by its binding ability in a functional ELISA.

Immobilized canine TNFRSF9-His at 10 µg/ml (100 µl/well) can bind rat Fc-TNFSF9 (Cat:503740), The EC₅₀ of rat Fc-TNFSF9 (Cat:503740) is 0.26-0.62 µg/ml.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ile 24

Molecular Mass

The recombinant canine TNFRSF9 comprises 173 amino acids and has a predicted molecular mass of 18.4 kDa. The apparent molecular mass of the protein is approximately 28-34 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Canine CD137 / 4-1BB / TNFRSF9 Protein (His Tag) SDS-PAGE