

# RPS13 antibody



Catalog Number: 114828

---

## Product name

RPS13 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

RPS13 Rabbit Polyclonal antibody. Positive WB detected in mouse uterus tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 14kd, 17kd

## Preparation

This antibody was obtained by immunization of RPS13 recombinant protein (Accession Number: NM\_001017). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse uterus tissue were subjected to SDS PAGE followed by western blot with Catalog No:114828(RPS13 antibody) at dilution of 1:500