

# SULT1B1 antibody



Catalog Number: 115828

## Product name

SULT1B1 antibody

## Specificity

human; other species not tested.

## Antibody description

SULT1B1 Rabbit Polyclonal antibody. Positive WB detected in human colon tissue, human liver tissue. Observed molecular weight by Western-blot: 35kd

## Preparation

This antibody was obtained by immunization of SULT1B1 recombinant protein (Accession Number: XM\_005265677). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:115828(SULT1B1 antibody) at dilution of 1:600