

Andecaliximab reference antibody



Catalog Number: 176828

Product name

Andecaliximab reference antibody

Specificity

Human MMP9

Antibody description

Anti-MMP9 Antibody (Andecaliximab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG4

Applications


ELISA

Validations




SDS-PAGE

Anti-MMP9 Antibody (andecaliximab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



SEC-HPLC

The purity of Anti-MMP9 Antibody (andecaliximab) is more than 95%, determined by



SEC-HPLC.

Bioactivity: ELISA

Immobilized human MMP9 His at 2 µg/mL can bind Anti-MMP9 Antibody (andecaliximab), EC50=0.006577 µg/mL.