

# CD3 delta antibody



Catalog Number: 107128

## Product name

CD3 delta antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CD3 delta Mouse Monoclonal antibody. Positive WB detected in Jurkat cells, A2780 cells. Observed molecular weight by Western-blot: 25-30 kDa

## Preparation

This antibody was obtained by immunization of CD3 delta recombinant protein (Accession Number: XM\_017018543). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107128(CD3D Antibody) at dilution of 1:500