

# Anti-Cdk4 (3F9) Mouse antibody



Catalog Number: 168128

## Product name

Anti-Cdk4 (3F9) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Cdk4 (3F9) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human CDK4 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB


## Dilutions

WB: 1/1000;Other 1/1000

## Validations



Western blot detection of CDK4 in HeLa,Lncap and F9 cell lysates using CDK4 mouse mAb (1:1000 diluted).Predicted band size:34KDa.Observed



band size:33KDa.

Western blot detection of CDK4 in HeLa,C6 and F9 cell lysates using CDK4 mouse mAb (1:1000 diluted).Predicted band size:34KDa.Observed band size:34KDa.