

# SPACA3 antibody



Catalog Number: 115528

## Product name

SPACA3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SPACA3 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse testis tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 28-32 kDa, 46-50 kDa

## Preparation

This antibody was obtained by immunization of SPACA3 recombinant protein (Accession Number: NM\_001317225). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115528(SPACA3 antibody) at dilution of 1:300



Immunofluorescent analysis of HEK-293 cells using Catalog No:115528(SPACA3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)