

UBE2J2 antibody



Catalog Number: 116528

Product name

UBE2J2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UBE2J2 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of UBE2J2 recombinant protein (Accession Number: XM_011540614). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116528(UBE2J2 Antibody) at dilution of 1:300