

MYH1 antibody

Catalog Number: 112928

Product name

MYH1 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

MYH1 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IP detected in mouse skeletal muscle tissue. Positive IHC detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 220-230 kDa

Preparation

This antibody was obtained by immunization of MYH1 recombinant protein (Accession Number: NM_005963). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, IP, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112928(MYH1 Antibody) at dilution of 1:1000



IP Result of anti-MYH1 (IP:Catalog No:112928, 5ug; Detection:Catalog No:112928 1:1000) with mouse skeletal muscle tissue lysate 4000ug.




Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:112928(MYH1 Antibody) at dilution of 1:50

MYH1 antibody

Catalog Number: 112928

(under 10x lens)



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:112928(MYH1 Antibody) at dilution of 1:50 (under 40x lens)