

Oxelumab (Synonyms: HUMAB OX 40L)



Catalog Number: 176928

Product name

Oxelumab (Synonyms: HUMAB OX 40L)

Specificity

Human TNFSF4 / OX40L / CD252

Antibody description

Anti-TNFSF4 / OX40L / CD252 Antibody (Oxelumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations



SDS-PAGE

Anti-TNFSF4 / OX40L / CD252 Antibody (oxelumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-TNFSF4 / OX40L / CD252 Antibody (oxelumab) is 99.04%, determined by SEC-HPLC.