

Human DAO/DAMOX (His tag) recombinant protein



Catalog Number: 519928

General Information

Gene Name Synonym

D-Amino-Acid Oxidase, DAAO, DAMOX, DAO

Protein Construction

A DNA sequence encoding the Human DAO/DAMOX Met1-Leu347 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM DTT, 2mM EDTA, pH 8.0 .

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human DAO/DAMOX (His tag) recombinant protein