

# CDK11B Polyclonal Antibody



Catalog Number: 162228

---

## Product name

CDK11B Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to CDK11B

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human CDK11B (NP\_277021.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


IHC, IF

## Dilutions

IHC 1:50 - 1:200

IF 1:50 - 1:200

## Validations



Immunofluorescence - CDK11B Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using CDK11B antibody at dilution of 1:100. Blue: DAPI for nuclear staining.