

# NQO2 antibody



Catalog Number: 113228

## Product name

NQO2 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

## Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

## Validations

### Antibody description

NQO2 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue. Positive IF detected in HeLa cells. Positive WB detected in HeLa cells, human liver tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 26kd



### Preparation

This antibody was obtained by immunization of NQO2 recombinant protein (Accession Number: NM\_001318940). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113228(NQO2 antibody) at dilution of 1:500

### Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

### Storage

Store at -20°C. DO NOT ALIQUOT

### Clonality

Polyclonal

### Ig Type

Rabbit IgG


### Applications

ELISA, WB, IHC, IF, IP


### Dilutions

Recommended Dilution:

WB: 1:1000-1:10000



Immunohistochemical of paraffin-embedded human brain using Catalog No:113228(NQO2 antibody) at dilution of 1:100 (under 10x lens)




# NQO2 antibody

Catalog Number: 113228

---

Immunohistochemical of paraffin-embedded human brain using Catalog No:113228(NQO2 antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HeLa cells, using NQO2 antibody Catalog No:113228 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). 

IP Result of anti-NQO2 (IP:Catalog No:113228, 3ug; Detection:Catalog No:113228 1:500) with HeLa cells lysate 1200ug.