

Anti-OR7C1 antibody



Catalog Number: 176228

Product name

Anti-OR7C1 antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IF

Antibody description

Rabbit polyclonal antibody to OR7C1

Dilutions

WB: 1:500-1:2000

Preparation

IF: 1:100-1:500

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or7c1. aa range:199-248

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Validations



Fig1:: Western Blot analysis of various cells using Olfactory receptor 7C1 Polyclonal Antibody diluted at 1: 500