

UPK1A antibody

Catalog Number: 116628

Product name

UPK1A antibody

Specificity

Human; other species not tested.

Antibody description

UPK1A Rabbit Polyclonal antibody. Positive IP detected in mouse bladder tissue. Positive IHC detected in human bladder tissue. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of UPK1A recombinant protein (Accession Number: NM_007000). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IP


Dilutions

Recommended Dilution:


IP: 1:200-1:1000

IHC: 1:20-1:200

Validations



IP Result of anti-UPK1A (IP:Catalog No:116628, 4ug; Detection:Catalog No:116628 1:300) with mouse bladder tissue lysate 3200ug.



Immunohistochemistry of paraffin-embedded human bladder tissue slide using Catalog No:116628(UPK1A Antibody) at dilution of 1:1600 (under 10x lens)