

Anti-LIFR/CD118 antibody



Catalog Number: 102328

Product name

Monoclonal

Anti-LIFR/CD118 antibody

Ig Type

Immunogen

Rabbit IgG

[Human LIFR/CD118 \(His Tag\) recombinant protein](#)

Applications

ELISA, IHC-P

Specificity

Dilutions

Human LIFR / CD118

ELISA: 0.1-0.2 µg/mL

Antibody description

Rabbit monoclonal to LIFR/CD118

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human LIFR. The detection limit for Human LIFR is 0.00975 ng/well.

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human LIFR / CD118 (rh LIFR; NP_001121143.1; Met 1-Ser 833).

IHC-P: 1-10 µg/mL

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Validations


Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality



LIFR / CD118 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human LIFR in human placenta with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections). Positive staining was localized to trophoblast.