

Anti-TAS2R14 antibody

Catalog Number: 176328

Product name

Anti-TAS2R14 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to TAS2R14

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tas2r14. aa range:229-278

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western blot analysis of lysates from MCF-7 cells, using TAS2R14 Antibody. The lane on the right is blocked with the synthesized peptide.




Fig2.; Immunofluorescence analysis of MCF7 cells, using TAS2R14 Antibody. The picture on the right is blocked with the synthesized peptide.