

H1N1 HA (His Tag) recombinant protein



Catalog Number: 503728

General Information

Protein Construction

A DNA sequence encoding the N-terminal segment (Met 1-Arg 343) of the influenza hemagglutinin (A/WSN/33 (H1N1)) (ACF54598.1), termed as HA1, was fused with a C-terminal polyhistidine tag.

Organism

H1N1

Expression Host

Human Cells

QC Testing

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 18

Molecular Mass

The secreted recombinant HA1 subunit of influenza A H1N1 HA (A/WSN/1933 (H1N1)) comprises 337 amino acids and has a predicted

molecular mass of 38 kDa. As a result of glycosylation, the apparent molecular mass of the protein is approximately 45-50 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H1N1 (A/WSN/1933) Hemagglutinin Protein (HA1 Subunit) (His Tag) SDS-PAGE