

RNASEH1 antibody

Catalog Number: 114728

Product name

RNASEH1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RNASEH1 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, human brain tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of RNASEH1 recombinant protein (Accession Number: NM_002936). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114728(RNASEH1 antibody) at dilution of



1:500

IP Result of anti-RNASEH1 (IP:Catalog No:114728, 3ug; Detection:Catalog No:114728 1:500) with HeLa cells lysate 1200ug.