

Anti-SF3B1 Rabbit antibody



Catalog Number: 167728

Product name

Anti-SF3B1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

SF3B1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human SF3B1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

FC: 1/200

Validations



Western blot detection of SF3B1 in HeLa, CHO-K1, C6, Ramos using SF3B1 Rabbit mAb(1:1000 diluted)