

GLT25D2 antibody



Catalog Number: 111038

Product name

GLT25D2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GLT25D2 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, Neuro-2a cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of GLT25D2 recombinant protein (Accession Number: NM_015101). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:111038(GLT25D2 Antibody) at dilution of 1:600