

# GLT25D2 antibody



Catalog Number: 111038

## Product name

Polyclonal

GLT25D2 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse; other species not tested.

## Applications

## Antibody description

ELISA, WB

GLT25D2 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, Neuro-2a cells. Observed molecular weight by Western-blot: 65 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of GLT25D2 recombinant protein (Accession Number: NM\_015101). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:111038(GLT25D2 Antibody) at dilution of 1:600