

Mouse Neuroserpin/Serpini1 (His Tag) recombinant protein



Catalog Number: 502038

General Information

Gene Name Synonym

Peptidase inhibitor 12; Serine protease inhibitor 17; Serpin I1

Protein Construction

A DNA sequence encoding the mouse SERPINI1 (O35684) (Met 1-Leu 410) was expressed, with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Thr 17

Molecular Mass

The secreted recombinant mouse SERPINI1

comprises 405 amino acids and has a calculated molecular mass of 46 kDa. The apparent molecular mass of rmSERPINI1 is approximately 45-50 kDa in SDS-PAGE under reducing conditions as a result of glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse SerpinI1 / Neuroserpin Protein (His Tag)
SDS-PAGE