

CD109 antibody



Catalog Number: 107038

Product name

CD109 antibody

Specificity

Human; other species not tested.

Antibody description

CD109 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 162 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: XM_005248659). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107038(CD109 Antibody) at dilution of 1:4000