

# Anti-Osteopontin Rabbit antibody



Catalog Number: 168038

## Product name

Anti-Osteopontin Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Osteopontin Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Osteopontin

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC

## Validations



Western blot detection of Osteopontin in HeLa,CHO-K1,C6,Jurkat cell lysates using Osteopontin Rabbit mAb(1:1000 diluted).Predicted band size:35kDa.Observed band size:35kDa.