

CCZ1 antibody

Catalog Number: 109038

Product name

CCZ1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CCZ1 Rabbit Polyclonal antibody. Positive IP detected in NIH/3T3 cells. Positive WB detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 50-56 kDa

Preparation

This antibody was obtained by immunization of CCZ1 recombinant protein (Accession Number: NM_015622). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



IP Result of anti-C7orf28A (IP:Catalog No:109038, 4ug; Detection:Catalog No:109038 1:500) with



NIH/3T3 cells lysate 1500ug.

NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109038(C7orf28A Antibody) at dilution of 1:1000