

RAB39B antibody

Catalog Number: 114438

Product name

RAB39B antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAB39B Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in human brain tissue, Y79 cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 25 kDa

Preparation

This antibody was obtained by immunization of RAB39B recombinant protein (Accession Number: NM_171998). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114438(RAB39B antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:114438(RAB39B antibody) at dilution of 1:50



(under 10x lens)

IP Result of anti-RAB39B (IP:Catalog No:114438,

RAB39B antibody



Catalog Number: 114438

4ug; Detection:Catalog No:114438 1:300) with mouse brain tissue lysate 4000ug.