

ORM2-Specific antibody



Catalog Number: 107438

Product name

ORM2-Specific antibody

Specificity

Human; other species not tested.

Antibody description

ORM2-Specific Mouse Monoclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in human plasma tissue, A549 cells. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of ORM2-Specific recombinant protein (Accession Number: NM_000608). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:107438(ORM2 Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol) fixed HepG2 cells using Catalog No:107438(ORM2 Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)