

Anti-TRIM33 (3H7) Mouse antibody



Catalog Number: 168438

Product name

Anti-TRIM33 (3H7) Mouse antibody

Specificity

Human

Antibody description

TRIM33 (3H7) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human TIF1 gamma / TRIM33 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2a


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of TIF1 gamma / TRIM33 in MCF7 and HeLa cell lysates using TIF1 gamma / TRIM33 mouse mAb (1:1000 diluted). Predicted band size: 120KDa. Observed band size: 140KDa.