

Human ERK2 / MAPK1 / MAPK2 (GST Tag) recombinant protein



Catalog Number: 500838

General Information

Gene Name Synonym

ERT1; Extracellular signal-regulated kinase 2; MAP kinase isoform p42; Mitogen-activated protein kinase 2

Protein Construction

A DNA sequence encoding the human ERK2 (NP_002736.3) (Met 1-Ser 360) was fused with the GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

No Kinase Activity

Purity

> 98 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human ERK2/GST chimera consists of 585 amino acids and predicts a molecular mass of 67 kDa as estimated in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 50mM Tris, 100mM NaCl, 0.5mM PMSF, 10% Glycerol, pH 8.0

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human ERK2 / MAPK1 / MAPK2 Protein (GST Tag) SDS-PAGE