

H7N9 Neuraminidase (Active) recombinant protein



Catalog Number: 504838

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/Shanghai/1/2013 (H7N9)) Neuraminidase (His36-Leu465) was expressed, the cell lysates are collected, and bio-activity was tested.

Organism

H7N9

Expression Host

Human Cells

QC Testing

Activity

Measured by its ability to cleave a fluorogenic substrate, 2'-(4-Methylumbelliferyl)- α -D-N-acetylneuraminic acid

The specific activity is > 20 U

One unit is defined as the amount of enzyme required to cleave 1 nmole of 2'-(4-Methylumbelliferyl)- α -D-N-acetylneuraminic acid per minute at pH 7.5 at 37°C.

Endotoxin

< 1.0 EU per μ g of the protein as determined by

the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Molecular Mass

The influenza H7N9 virus Neuraminidase comprises 448 amino acids

Formulation

Lyophilized from sterile PBS, 0.6% Triton-100, 6% trehalose 5.3% mannitol PH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

It is recommended that 1 ml sterile water be added to the vial to prepare a stock solution.