

AGTR1 antibody

Catalog Number: 107838

Product name

AGTR1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AGTR1 Rabbit Polyclonal antibody. Positive WB detected in PC-12 cells, mouse kidney tissue. Positive IHC detected in human heart tissue, human liver cancer tissue. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of AGTR1 recombinant protein (Accession Number: NM_004835). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



PC-12 cells were subjected to SDS PAGE followed by western blot with Catalog No:107838(AGTR1 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:107838(AGTR1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:107838(AGTR1 Antibody) at dilution of 1:50

AGTR1 antibody



Catalog Number: 107838

(under 40x lens)