

# Human GGPS1 (His tag) recombinant protein



Catalog Number: 518838

## General Information

### Gene Name Synonym

Geranylgeranyl Pyrophosphate Synthase, GGPP Synthase, GGPPSase, (2E,6E)-Farnesyl Diphosphate Synthase, Dimethylallyltranstransferase, Farnesyl Diphosphate Synthase, Farnesyltranstransferase, Geranylgeranyl Diphosphate Synthase, Geranyltranstransferase, GGPS1

### Protein Construction

A DNA sequence encoding the Human GGPS1 Met1-Glu300 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0 .

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human GGPS1 (His tag) recombinant protein