

SF3B3 antibody

Catalog Number: 115138

Product name

SF3B3 antibody

WB: 1:1000-1:10000

IP: 1:2000-1:20000

Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:20-1:200

Antibody description

SF3B3 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IHC detected in human gliomas tissue. Positive IP detected in rat brain tissue. Positive WB detected in HeLa cells, human brain tissue, mouse heart tissue, rat brain tissue, rat heart tissue. Observed molecular weight by Western-blot: 130kd

Validations



Preparation

This antibody was obtained by immunization of SF3B3 recombinant protein (Accession Number: BC003146). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115138(SF3B3 antibody) at dilution of 1:10000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IF, IP

Dilutions


Recommended Dilution:




Immunohistochemical of paraffin-embedded human gliomas using Catalog No:115138(SF3B3 antibody) at dilution of 1:100 (under 40x lens)

SF3B3 antibody

Catalog Number: 115138



Immunofluorescent analysis of HeLa cells, using SF3B3 antibody Catalog No:115138 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent



DNA dye).

IP Result of anti-SF3B3 (IP:Catalog No:115138, 4ug; Detection:Catalog No:115138 1:5000) with rat brain tissue lysate 4800ug.