

Anti-TLR6 antibody



Catalog Number: 106138

Product name

Anti-TLR6 antibody

Specificity

Human TLR6

Antibody description

Rabbit polyclonal to TLR6

Preparation

Produced in rabbits immunized with A synthetic peptide corresponding to the center region of the human CD286, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 0.1-2 µg/mL

Validations



CD286 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human TLR6 in human hepatoma with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded



sections).

CD286 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human TLR6 in human breast carcinoma with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded sections).