

# Melrilimab reference antibody



Catalog Number: 177138

## Product name

Melrilimab reference antibody

## Specificity

Human IL-1RL1 / ST2 / IL-33R

## Antibody description

Anti-IL-1RL1 / ST2 / IL-33R Antibody (Melrilimab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG2SA

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-IL-1RL1 / ST2 / IL-33R Antibody (melrilimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-IL-1RL1 / ST2 / IL-33R Antibody (melrilimab) is more than 95%, determined by SEC-HPLC.