

Anti-GAB1 Mouse antibody



Catalog Number: 169138

Product name

Anti-GAB1 Mouse antibody

Specificity

Human

Antibody description

GAB1 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human GAB1 expressed in E. Coli. ;

Formulation

Purified antibody in PBS with 0.05% sodium azide and 0.5% protein stabilizer ;

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

WB, IHC-P, FC

Validations




Immunohistochemical analysis of paraffin-embedded kidney tissues using GAB1 mouse mAb with DAB staining.



with DAB staining.



Western blot analysis using GAB1 mouse mAb against HEK293 (1) cell lysate.



Flow cytometric analysis of Jurkat cells using GAB1 mouse mAb (green) and negative control (red).

Anti-GAB1 Mouse antibody

Catalog Number: 169138



Western blot analysis using GAB1 mAb against human GAB1 (AA (red).

Immunohistochemical analysis of paraffin-embedded bladder cancer tissues using GAB1 mouse mAb with DAB staining.