

SPAM1 antibody



Catalog Number: 115538

Product name

SPAM1 antibody

Specificity

Human; other species not tested.

Antibody description

SPAM1 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, human testis tissue, MCF7 cells. Observed molecular weight by Western-blot: 45-47 kDa

Preparation

This antibody was obtained by immunization of SPAM1 recombinant protein (Accession Number: NM_003117). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115538(SPAM1 antibody) at dilution of 1:500