

Human CA7 (His tag) recombinant protein



Catalog Number: 518538

General Information

Gene Name Synonym

Carbonic Anhydrase 7, Carbonate Dehydratase VII, Carbonic Anhydrase VII, CA-VII, CA7

Protein Construction

A DNA sequence encoding the Human CA7 Met1-Ala264 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0 .

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human CA7 (His tag) recombinant protein