

Human JNK1 / MAPK8 (GST Tag) recombinant protein



Catalog Number: 502938

General Information

Protein Construction

A DNA sequence encoding the full length of human MAPK8 isoform JNK1 alpha 2 (NP_620637.1) (Met 1-Arg 427) was expressed with the GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

No Kinase Activity

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human MAPK8/GST chimera

consists of 652 amino acids and predicts a molecular mass of 75 kDa. It migrates as an approximately 65 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 50mM Tris, 100mM NaCl, pH 8.0, 25% glycerol

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human JNK1 / MAPK8 Protein (GST Tag) SDS-PAGE