

RUNX2 antibody



Catalog Number: 114938

Product name

RUNX2 antibody

Specificity

Human, Rat; other species not tested.

Antibody description

RUNX2 Rabbit Polyclonal antibody. Positive WB detected in ROS1728 cells. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004348).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



ROS1728 cells were subjected to SDS PAGE followed by western blot with Catalog No:114938(RUNX2 Antibody) at dilution of 1:600