

Anti-CYP2A7 antibody

Catalog Number: 175938

Product name

Anti-CYP2A7 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CYP2A7

Preparation

This antigen of this antibody was synthesized peptide derived from the c-terminal region of human cyp2a7.

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.: Western Blot analysis of various cells using CYP2A7 Polyclonal Antibody diluted at 1: 2000