

NLE1 antibody



Catalog Number: 113238

Product name

NLE1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NLE1 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells. Observed molecular weight by Western-blot: 45-55 kDa

Preparation

This antibody was obtained by immunization of NLE1 recombinant protein (Accession Number: NM_018096). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113238(NLE1 Antibody) at dilution of 1:300