

PRPS1L1 antibody



Catalog Number: 114238

Product name

PRPS1L1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PRPS1L1 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 35kd

Preparation

This antibody was obtained by immunization of PRPS1L1 recombinant protein (Accession Number: NM_175886). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114238(PRPS1L1 antibody) at dilution of 1:1000