

# ADI-65534 antibody



Catalog Number: 177238

## Product name

ADI-65534 antibody

## Antibody description

ADI 65534 an neutralizing antibody against ANDV (New World Hantavirus) and PUUV (Old World Hantavirus)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

Lyophilized, PBS(pH7.2)

## Storage

Store at -20°C. Avoid freeze / thaw cycles.


## Clonality

Monoclonal

## Ig Type


Human IgG

## Validations



## SDS-PAGE

Lane 1: Reduced sample, Lane 2: Non-reduced sample



## SEC-HPLC

HPLC Purity: 98.75%