

ADI-65534 antibody



Catalog Number: 177238

Product name

ADI-65534 antibody

Antibody description

ADI 65534 an neutralizing antibody against ANDV (New World Hantavirus) and PUUV (Old World Hantavirus)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

Lyophilized, PBS(pH7.2)

Storage

Store at -20°C. Avoid freeze / thaw cycles.


Clonality

Monoclonal

Ig Type


Human IgG

Validations



SDS-PAGE

Lane 1: Reduced sample, Lane 2: Non-reduced sample



SEC-HPLC

HPLC Purity: 98.75%